Unique design

- To optimize the comfort and durability of our BraceID[®] walkers we've added a unique, new pump system to the air liner.
- The concept consists of a fully integrated pump and valve system that is directly connected to the foam liner.
- We are using flexible and durable materials which result in a system that keeps working in every condition.
- The new pump handles easily and is easy to lift up from, or attach it back onto the top strap.



Integrated Value system

Flexible and

Durable material

Easy to operate



-2

Features and Benefits

- Rigid and lightweight Walker with front and back panel for more complete protection and support.
- Rocker bottom design promotes natural gait.
- Inflatable AIR liner for better fixation of foot and ankle with increased comfort. Intermittent pneumatic compression not only promotes callus formation over fracture site but also reduces local swelling.

Indications Low Version

- Acute ankle sprain and soft tissue injuries of the lower leg
- Stress fractures and stable or internally fixed fractures of the foot and ankle
- Achilles tendon repair

Size	Shoe size	Standard model	Air model
Small	36 - 39	12-0001-S	12-0009-S
Medium	39 - 43	12-0001-M	12-0009-M
Large	43 - 46	12-0001-L	12-0009-L
X-Large	+ 46	12-0001-XL	12-0009-XL



Features and Benefits

- Supplied with additional pieces of foam padding to increase patient comfort.
- Heel wedges ⁽¹⁾ applied inside the Walker allow to position the ankle in plantar flexion (5 pieces).

Indications Low Version

- Immobilisation after foot or toe chirgurie
- Severe ankle sprains
- Metatarsal fractures
- Front or mid foot injuries



LOW Shield Walker

	Size	Shoe size	Standard model	Air model	
	Small	36 - 39	12-0008-S	12-0010-S	
	Medium	39 - 43	12-0008-M	12-0010-M	
	Large	43 - 46	12-0008-L	12-0010-L	
	X-Large	+ 46	12-0008-XL	12-0010-XL	

Product Composition Polyester, PU foam, Nylon, Fibre Glass, Plastic

