

ALFA-Flex[®] foot abduction brace



The ALFA-Flex foot abduction brace is used for the treatment of idiopathic clubfoot according to Dr Ponseti's conservative therapy approach. It is used in conjunction with a foot support as an orthosis to hold the corrected foot in the desired position after plaster treatment. Depending on the doctor's recommendation, this can last another three to four years - thus preventing a recurrence.

Abduction and dorsiflexion - all angles accurate to the degree

The ALFA-Flex foot abduction brace enables the precise adjustment of the angles required for Ponseti treatment.

Dorsiflexion between 0° and 20°

Abduction between 0° and 80°

Grows with the patient - one brace for the entire treatment

The ALFA-Flex foot abduction brace is precisely adjustable in length thanks to a slotted hole system up to a width of 17,5 cm to 30 cm. It can be used for the entire duration of therapy.

ALFA Flex XL Brace and ALFA-Flex XL centre oval

Brace and accessories are available for particularly tall children, with a width of 18.5 cm to 37 cm.

The matching foot supports

The foot supports (available separately) are connected to the foot abduction system with a handy click system. The ALFA-Flex brace is suitable for use with all Sameda shoe models.

Information and ordering

ALFA-Flex Brace, article no.: 0700-000

Order online at www.semeda-shop.de

WE 5002-001